

KC4 Delta and KC1 Delta

Semi Automated Coagulation Analyser

KC4 Delta and KC1 Delta are semi automated coagulation analysers with four or one test position(s), respectively, providing operators with a compact easy to use system. KC Delta series instruments use micro-mechanical clot detection technology for clotting assays.



Technology

- “Gold standard” mechanical detection
- Pipette auto start testing
- LCD display and optional printing of results
- Programmable test modes, single or duplicate testing



Measuring Features

- Store reagent ISI values for automatic INR calculation including calibration curves
- Preparation and incubation area for samples and reagents
- Suitable for STAT and routine testing
- Test menu for PT, APTT, Fibrinogen, Factors
- Maintenance free operation

